

Oratect®

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
highly specific enzymatic reaction:
Alcohol Oxidase H

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Fig. f Example of presumptive positive alcohol test

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During testing, oral fluid is collected on the alcohol pad and saturates the alcohol pad. If no alcohol is present in the oral fluid, the alcohol pad remains colorless (remains white or cream color) because there is no alcohol in the oral fluid to react with enzymes to start the color reaction. If alcohol is present in the oral fluid, the alcohol pad changes to green or blue color because the alcohol reacts with alcohol oxidase to produce aldehyde and peroxide. The peroxide reacts with peroxidase in the presence of hydrogen donor to produce a blue color. Therefore, the presence of green to blue color at the alcohol pad window indicates a presumptive positive result for alcohol.

(2) Alcohol Test: The alcohol pad contains Tetramethylbenzidine, Alcohol Oxidase, Peroxidase, Buffer and Stabilizing Proteins

Fig. e Example of negative alcohol test

Precision – Alcohol Test For the alcohol test, saliva was obtained by rinsing with positive ethanol control solutions at various B.A.C. (0.02%, 0.04% and 0.06%). Negative saliva was used to test at 0% concentration. For each concentration, a total of 15 tests were performed to validate the test performance. The results of the OratectPlus™ Oral Fluid Drug and Alcohol Screen Device are summarized below:

Test	Total # of Test/ Concentration	B.A.C.							
		0.00%		0.02%		0.04%		0.06%	
		-	+	-	+	-	+	-	+
Alcohol	15	15	0	1	14	0	15	0	15

English

(2) Alcohol Test: A pad coated with enzymes, turns to color shades of green and blue on contact with alcohol in the oral fluids. The alcohol pad employs a solid phase chemistry which uses the following

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Produttore: Branan Medical Corp. – Irvine (CA) USA
No. Reg. 3006026304 _F.D.A.(U.S.A.)